




105-RS1-Rack




105-RS3-Draw

## SU105UE - Rack Options

105-RS1-Rack and  
105-RS3-Draw



105-RS2-Rack and  
105-RS4-Draw



### Complete Rack Solutions – Includes racks, boxes, and dividers

Part Number	Box Size	Box Capacity	Rack	Layout
105-RS1-Rack	2"	72	Side access rack system to complete the SU105UE undercounter freezer. 2" boxes (6 racks, 72 box capacity), <b>BOXES &amp; DIVIDERS INCLUDED</b>	2 compartments, 1 shelf 3 racks per shelf
105-RS2-Rack	3"	54	Side access rack system to complete the SU105UE undercounter freezer. 3" boxes (6 racks, 54 box capacity), <b>BOXES &amp; DIVIDERS INCLUDED</b>	2 compartments, 1 shelf 3 racks per shelf
105-RS3-Draw	2"	72	Drawer style rack system to complete the SU105UE undercounter freezer. 2" boxes (6 racks, 72 box capacity), <b>BOXES &amp; DIVIDERS INCLUDED</b>	2 compartments, 1 shelf 3 racks per shelf
105-RS4-Draw	3"	54	Drawer style rack system to complete the SU105UE undercounter freezer. 3" boxes (6 racks, 54 box capacity), <b>BOXES &amp; DIVIDERS INCLUDED</b>	2 compartments, 1 shelf 3 racks per shelf

### Individual Side Access Racks

Part Number	Box Size	Description
105-R-342	2"	Side access rack for the SU105UE undercounter freezer. 2" boxes (12 box capacity, 4 high x 3 deep), <b>RACK ONLY</b>
105-R-342-81	2"	Side access rack for the SU105UE undercounter freezer. 2" boxes (12 box capacity, 4 high x 3 deep), <b>BOXES &amp; DIVIDERS INCLUDED</b>
105-R-333	3"	Side access rack for the SU105UE undercounter freezer. 3" boxes (9 box capacity, 3 high x 3 deep), <b>RACK ONLY</b>
105-R-333-81	3"	Side access rack for the SU105UE undercounter freezer. 3" boxes (9 box capacity, 3 high x 3 deep), <b>BOXES &amp; DIVIDERS INCLUDED</b>

### Individual Drawer Racks

Part Number	Box Size	Description
105-DR-342	2"	Drawer style rack for the SU105UE undercounter freezer. 2" boxes (12 box capacity, 4 high x 3 deep), <b>RACK ONLY</b>
105-DR-342-81	2"	Drawer style rack for the SU105UE undercounter freezer. 2" boxes (12 box capacity, 4 high x 3 deep), <b>BOXES &amp; DIVIDERS INCLUDED</b>
105-DR-333	3"	Drawer style rack for the SU105UE undercounter freezer. 3" boxes (9 box capacity, 3 high x 3 deep), <b>RACK ONLY</b>
105-DR-333-81	3"	Drawer style rack for the SU105UE undercounter freezer. 3" boxes (9 box capacity, 3 high x 3 deep), <b>BOXES &amp; DIVIDERS INCLUDED</b>

### Boxes & Dividers

Part Number	Box Size	Description
105-CB2-81	2"	2" Cardboard Box with 81 place divider, Box dimensions are 5.25"W x 5.25"D x 2"H
105-CB3-81	3"	3" Cardboard Box with 81 place divider, Box dimensions are 5.25"W x 5.25"D x 3"H
D81	2" or 3"	81 Place Divider for 2" or 3" boxes