





FLUIDX SAMPLE STORAGE SOLUTIONS



STARTER PACKS - SAVE UP TO 20% OFF

CATALOGUE NUMBER	DESCRIPTION	PRICE (ex GST)
FX-LITE	FluidX LITE Package includes: 1 x Ziath Uno Single Tube Readers (ZTS-UNO-80) 1 x EzyCap Decapper Recapper (EzyCap) 1 x Manual Single Tube Picker of your choice (10-5010, 10-5020) 1 x Case of 2D Coded Tubes of your choice	\$3,685 
FX-PLUS	FluidX PLUS Package includes: 1 x Ziath Uno Single Tube Readers (ZTS-UNO-80) 1 x EzyCap Decapper Recapper (EzyCap) 1 x IntelliXcap™ Semi-Automated Handheld Screw Cap Decapper 1 x Manual Single Tube Picker of your choice (10-5010, 10-5020) 1 x Case of 2D Coded Tubes of your choice	\$7,744 
FX-STARTER	FluidX STARTER Package includes: 1 x Ziath Uno Single Tube Readers (ZTS-UNO-80) 1 x Ziath DP5 Mirage Rack Reader (DP5-M-STD-80) 1 x IntelliXcap™ Semi-Automated Handheld Screw Cap Decapper 1 x Manual Single Tube Picker of your choice (10-5010, 10-5020) 1 x Case of 2D Coded Tubes of your choice	\$17,626 
FX-PREMIUM	FluidX PREMIUM Package includes: 1 x Ziath Uno Single Tube Readers (ZTS-UNO-80) 1 x Ziath DP5 Mirage Rack Reader (DP5-M-STD-80) 1 x Semi-Automated Screw Cap Decapper/Recapper, (Aperio™) 1 x Manual Single Tube Picker of your choice (10-5010, 10-5020) 1 x Case of 2D Coded Tubes of your choice	\$25,294 

Valid to 31 December 2024. Subject to Terms and Conditions.



ZTS-UNO-80



DP5-M-STD-80



IntelliXcap™ Semi-Automated



Aperio™



10-5010, 10-5020



EzyCap

Anatomy of an FluidX Tube External Thread, Tri-Coded

The Tri-Coded tube with external thread cap offers our most secure seal. This highly secure design also offers significant benefits over internal thread caps including, a higher working volume.

External Threaded Cap:

- Automation friendly
- Our most secure cap
- Enables greater working volume
- Designed to prevent over tightening

High Quality Virgin Polypropylene:

- No detectable leachables or extractables
- Manufactured from USP Class VI resin in a Class 8 clean room
- Endotoxin, DNA/RNA, heavy metals & animal free

Side Coding:

- Linear Barcode
- Human Readable Number
- Permanent laser etching

Clear Window:

- Allows visual inspection
- Optional space to print direct due to high quality manufacturing

2D Coded:

- Enables whole rack or single tube reading
- High contrast enabling reliable reading
- Permanent laser etching

