











# CoolCell® Containers

## Controlled Rate Alcohol-Free Cell Freezing



- ✓ Consistent freeze rate
- ✓ No ongoing costs
- ✓ Safe and easy to use



PROMOTIONAL OFFER PRODUCTS	PRICE (EX GST)	ALCOHOL-FILLED CONTAINER	ALCOHOL-FREE CELL FREEZING CONTAINER
<b>CoolCell LX</b> Cell freezing container for 12 x 1ml to 2ml cryo tubes • BCS-405  • BCS-405O  • BCS-405PK 	<b>\$315</b>	<b>Requires isopropanol</b> <ul style="list-style-type: none"> <li>• Replace alcohol every 5 uses</li> <li>• Keep track of number of uses</li> <li>• Pre-cool alcohol in refrigerator</li> </ul>	<b>No alcohol</b> <ul style="list-style-type: none"> <li>• No fluids</li> <li>• No pre-cooling</li> </ul>
<b>CoolCell FTS30</b> Cell freezing container for 30 x 1ml to 2ml cryo tubes • BCS-170  • BCS-170G  • BCS-170O  • BCS-170PK 	<b>\$860</b>	<b>Inconsistent freeze rate</b> <ul style="list-style-type: none"> <li>• Alcohol degradation induces variability</li> <li>• Two circles of wells - different freeze rates for each</li> </ul>	<b>No variability</b> <ul style="list-style-type: none"> <li>• All vials have uniform freeze rate</li> <li>• Radially symmetric design ensures consistency</li> </ul>
<b>CoolCell® 1ml FX</b> Cell freezing container for 12 x 1ml 96-format Sample Tubes • BCS-407P  • BCS-407PO 	<b>\$315</b>	<b>Approximately \$350/year</b> <ul style="list-style-type: none"> <li>• Change alcohol weekly</li> <li>• Disposal of hazardous waste</li> </ul>	<b>No on-going cost</b> <ul style="list-style-type: none"> <li>• No alcohol purchase or disposal</li> </ul>
<b>CoolCell 5ml LX</b> Cell freezing container for 12 x 3.5ml to 5ml cryo tubes • BCS-406 	<b>\$420</b>	<b>Difficult to handle</b> <ul style="list-style-type: none"> <li>• Screw cap difficult to remove when frozen</li> <li>• Frozen unit is slippery and cold to touch</li> </ul>	<b>No stuck lids</b> <ul style="list-style-type: none"> <li>• Lid comes off easily when frozen</li> <li>• Not cold to the touch when removing from the -80°C freezer</li> </ul>
		<b>Wait between runs</b> <ul style="list-style-type: none"> <li>• Takes &gt;1 hour for device to return to room temperature for re-use</li> </ul>	<b>Quick re-use time</b> <ul style="list-style-type: none"> <li>• Ready to use again after five minutes</li> </ul>
		<b>Large thermal mass impacts local freezer area</b> <ul style="list-style-type: none"> <li>• Large heat capacity removed from alcohol impacts nearby samples</li> </ul>	<b>Low impact on freezer</b> <ul style="list-style-type: none"> <li>• 1/3 the heat impact on freezer compared to alcohol-filled units</li> </ul>

# TruCool® Accessories


## TruCool® Ice Pans and Buckets

- ✓ Non-toxic
- ✓ Recyclable ethyl-vinyl acetate (EVA)
- ✓ Suitable for use with ice, dry ice, liquid nitrogen, alcohol slurries



### PROMOTIONAL PRICING

Size	Dimensions (L x W x H)	Purple	Lime green	Orange	Pink	Blue	Green	Red	PRICE (ex GST)
Mini 1L Ice Pan	15.9 x 15.9 x 9.5 cm	BCS-211PL	BCS-211GR	BCS-211OR	BCS-211PK	BCS-211B	BCS-211	BCS-212	<b>\$120</b>
Midi 4L Ice Pan	31.1 x 22.2 x 11.4 cm	BCS-113PL	BCS-113GR	BCS-113OR	BCS-113PK	BCS-113B	BCS-113	BCS-114	<b>\$194</b>
Midi 4L Ice Pan with Lid	35.0 x 22.5 x 11.4cm	BCS-117PL	BCS-117GR	BCS-117OR	BCS-117PK	BCS-117B			<b>\$224</b>
Maxi 9L Ice Pan	40.6 x 31.8 x 11.4 cm	BCS-111PL	BCS-111GR	BCS-111OR	BCS111PK	BCS-111B	BCS-111	BCS-112	<b>\$259</b>
Maxi 9L Ice Pan with Lid	44.5 x 32.5 x 11.4 cm	BCS-118PL	BCS-118GR			BCS-118B			<b>\$300</b>
Round 2.5L Ice Bucket with Lid	Top diameter: 24.1 cm, Height: 12.1 cm	BCS-115-25PL	BCS-115-25GR	BCS-115-25OR	BCS-115-25PK	BCS-115-25B	BCS-115-25G	BCS-115-25R	<b>\$148</b>
Round 4L Ice Bucket with Lid	Top diameter: 26.0 cm, Height: 18.4 cm	BCS-115PL	BCS-115GR	BCS-115OR	BCS115PK	BCS-115B	BCS-115	BCS-115R	<b>\$164</b>

	FREEZER STORAGE				CRYOSTORAGE			
	2-inch Hinged CryoBox with adjustable plastic grid, 81-place		2-inch Hinged CryoBox with fixed plastic insert, 100-place		3.5-inch Hinged CryoBox with adjustable plastic grid, 81-place		2-inch Hinged CryoBox with adjustable plastic grid, LN2 drain holes, 81-place	
<b>Qty per pack</b>	5	50	5	50	6	30	5	50
<b>PRICE (EX GST)</b>	<b>\$70</b>	<b>\$580</b>	<b>\$75</b>	<b>\$546</b>	<b>\$90</b>	<b>\$408</b>	<b>\$63</b>	<b>\$546</b>
<b>Green</b>	BCS-206G	BCS-207G	BCS-209G	BCS-220G	BCS-215G	BCS-219G	BCS-217G	BCS-221G
<b>Purple</b>	BCS-206P	BCS-207P	BCS-209P	BCS-220P	BCS-215P	BCS-219P	BCS-217P	BCS-221P
<b>Multi-Colour</b>	BCS-206MC							
<b>Dimensions (L x W x H)</b>	13.3 x 13.3 x 5.1 cm 5.25 x 5.25 x 2.0 in		13.3 x 13.3 x 5.1 cm 5.25 x 5.25 x 2.0 in		13.3 x 13.3 x 9.2 cm 5.25 x 5.25 x 3.5 in		13.3 x 13.3 x 5.1 cm 5.25 x 5.25 x 2.0 in	

Valid to 31<sup>st</sup> December 2024. Subject to Terms and Conditions