

Select Your Bullet Blender

Start with sample size

Choose the Bullet Blender that's right for your lab



up to 0.3g or 0.3 mL
1.5 mL tubes
24 samples



up to 1 g or 1 mL
5 mL tubes
12 samples



up to 3 g or 3 mL
50 mL tubes
8 samples



Storm Pro
small footprint



5E Pro
small footprint



Gold
Maintains samples at 4°C
Very quiet operation



5E Gold
Maintains samples at 4°C
Very quiet operation



50 Gold
Maintains samples at 4°C
Very quiet operation



Have any questions or need more help?

Other models are available for special purposes.
Contact your sales representative for further assistance.