



# EzyCap™

*Electronic Decapper/Recapper and Tube Picker  
A BioTools Exclusive*

## EzyCap™

Developed by BioTools – The EzyCap™ is a lightweight, versatile and easy to use tool specifically designed to decap and recap FluidX 2D barcode tubes and enable the user to individually pick and place tubes from within 49/96 way racks and 10 x 10 Cryoboxes.

## Key Features

- » Fully Electronic Operation
- » Minimises Risk of RSI
- » Delivers consistent sealing torque
- » Supports Range of Tube Sizes: 0.3, 0.7, 1.0 and 2.0mL External Threaded Screw Cap Tubes
- » Tube Picking – *Enables user to pick and place individual tubes*
- » Lightweight and Ergonomic

## EzyCap™ Kit Includes

- 1 EzyCap™ Driver
- 2 0.3, 0.7, 1.0 ml Tube Adapter
- 3 2.0 ml Tube Adapter
- 4 AAA Batteries
- 5 Mounting Bracket
- 6 Protective Carry Case

