







# IntelliXcap Automated Septum Cap Decapper/Recapper Adapter Matrix

The IntelliXcap™ Septum Cap Decapper/Recapper includes three adaptors and locking teeth for use with Azenta sample tubes. Extra adaptors can be purchased separately, as well as adaptors for use with other manufacturer's tubes.

The matrix table below shows which adaptors to use with each tube type.

Rack Part Number	Rack Description	Caps Part Number	Rack Adapter	Rack Adapter Spacer	Locking Teeth	Locking Teeth Spacers	
<b>Azenta Sample Tubes</b>							
<i>Adapters included with IntelliXcap Septum Caps (46-5030)</i>							
66-62321-L	96-format 0.9 ml, Internal Thread, Next-Gen Dual-coded (uncapped, high base rack with Tube Lock)	65-73000	A	B	1	2	
66-62325	96-format 0.3 ml, Internal Thread, Next-Gen Dual-coded (uncapped, low base rack)	65-73000	A	D	1	3	
67-0753-01	96-format 0.48 ml, Internal Thread, Next-Gen Jacket Tri-coded (uncapped, low base rack)	65-73000	A	C	1	2	
<b>Other Manufacturer's Tubes</b>							
<i>Adapters available separately</i>							
3791	1.4 ml Matrix 2D Tubes	4465	A	E	1	4	 46-5051