



Manufacturer	Alpaqua Engineering, LLC	
Part #	A000350	
Manufacturer's Website	http://www.alpaqua.com/Products/MagnetPlates/LEMagnetPlate.aspx	
Description	Plastic/Metal Plate with 96 Ring Magnets for Magnetic Bead Separations	
Country of Origin	U.S.	
Harmonized Code	8505110070 - Sintered neodymium-iron-boron, permanent magnets & articles intended to become permanent magnets after magnetization	
US Export Restrictions	ECCN for this item is 0A999 (Ring Magnets). Generally, no export license from the U.S. Department of Commerce - Bureau of Industry and Security (BIS) is required except for export to North Korea (License requirement: AT (Anti-Terrorism)). Other restrictions governing export or re-export to denied countries, entities, and individuals may apply. See http://export.gov for details.	
Materials	NdFeB Magnets, Acetal Copolymer, Aluminum, Stainless Steel	
Dimensions	W 127.8 mm x L 85.6 mm x H 21.2 mm	
Weight	Net Weight:	360 g = 0.79 lb
	Shipping Weight:	400 g = 0.88 lb
MSDS	Available upon request	
RoHS/REACH	This product is compliant with RoHS/REACH regulations in effect on the date shown below.	
Shipping Restrictions	One single magnet plate in standard packaging has been tested for compliance with IATA regulations regarding transport by air of magnetized material. At 7 ft. of distance in any direction from the package the magnetic field was found to be safely below the threshold strength of 2 mGauss. One magnet plate is not considered a Dangerous Good and can be shipped by air without special labeling. (Test should be repeated for packages containing multiple magnet plates.)	
Additional Comments	Maximum Temperature: 80°C	

A000350 Product Passport.Docx - Last update: May 9, 2018