



Erlenmeyer Flask, 125ml, PETG

Code 020607-0125

Erlenmeyer Flask, 125ml, PETG

The 125ml Erlenmeyer Flask is ideal for the suspension of cell culture, media preparation, mixing and storage.

Manufactured for virgin PETG, the flask is robust and break resistant, thus increasing personal safety. The single use, flat bottomed flask is supplied with a vented HDPE screw cap and a 0.22µm microporous hydrophobic PTFE membrane for continuous gas exchange without contamination.

The flask and cap are packaged individually in a sterile zip sealable bag.



Product Features

- Erlenmeyer Flask, 125ml, PETG
- HDPE Screw Cap, Vented, Orange
- Carton of 50, Individually Bagged
- Gamma Radiation Sterilised
- Laboratory Use Only
- Disposable

Product Specifications

Flask Material	Virgin PETG
Flask Volume	125ml
Flask Bottom Type	Flat
Cap Type	Vent
Cap Material	Virgin HDPE
Cap Colour	Orange
Sealing Type	Screw Cap
Hydrophobic Membrane	0.22µm
Sterility	Gamma Radiation
Disposable	Yes
Certifications	ISO 13485 ISO 11137 ISO 90001

PRODUCT CODE

020607-0125



Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com

Packaging details

Pieces Per Bag	1
Bags Per Carton	50
Carton Dimensions	350mm x 420mm x 200mm
Carton Weight	3kg

Ordering Information

Code	Description
020607-0125	Erlenmeyer Flask, 125ml, PET, Orange Vent Cap, Sterile, 1/Bag, Carton of 50



Australia
1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada
1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand
+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States
646 568 5391
sales.us@westlab.com
www.westlab.com