








TruCool® Accessories



TruCool® Ice Pans and Buckets

Non-toxic, recyclable ethyl-vinyl acetate (EVA) foam containers for use with ice, dry ice, liquid nitrogen, alcohol slurries. Will not sweat, leak or skid on bench.



Size	Dimensions (L x W x H)	● purple	● lime green	● orange	● pink	● blue	● green	● red
 Mini 1L Ice Pan	15.9 x 15.9 x 9.5 cm	BCS-211PL	BCS-211GR	BCS-211OR	BCS-211PK	BCS-211B	BCS-211	BCS-212
 Midi 4L Ice Pan	31.1 x 22.2 x 11.4 cm	BCS-113PL	BCS-113GR	BCS-113OR	BCS-113PK	BCS-113B	BCS-113	BCS-114
 Midi 4L Ice Pan with Lid	35.0 x 22.5 x 11.4cm	BCS-117PL	BCS-117GR	BCS-117OR	BCS-117PK	BCS-117B		
 Maxi 9L Ice Pan	40.6 x 31.8 x 11.4 cm	BCS-111PL	BCS-111GR	BCS-111OR	BCS-111PK	BCS-111B	BCS-111	BCS-112
 Maxi 9L Ice Pan with Lid	44.5 x 32.5 x 11.4 cm	BCS-118PL	BCS-118GR			BCS-118B		
 Round 2.5L Ice Bucket with Lid	Top diameter: 24.1 cm, Height: 12.1 cm	BCS-115-25PL	BCS-115-25GR	BCS-115-25OR	BCS-115-25PK	BCS-115-25B	BCS-115-25G	BCS-115-25R
 Round 4L Ice Bucket with Lid	Top diameter: 26.0 cm, Height: 18.4 cm	BCS-115PL	BCS-115GR	BCS-115OR	BCS-115PK	BCS-115B	BCS-115	BCS-115R

TruCool® Hinged CryoBoxes

Patented hinged lid offers convenience and archival integrity, ensuring markings and vials remain in sync. Lid stays attached to base minimizing risk of separation and lid contamination. Lid is easy to open when frozen. Available in 9x9, 10x10, and vapor phase LN₂ compatible formats. Plastic 81-place grid has adjustable slats to accommodate multiple vial types. 2-inch box holds 1.0 mL or 2.0 mL cryogenic vials and microfuge tubes. 3.5-inch box holds 3.0 mL to 5.0 mL cryogenic vials.



	Freezer Storage				CryoStorage			
	2-inch Hinged CryoBox with adjustable plastic grid, 81-place		2-inch Hinged CryoBox with fixed plastic insert, 100-place		3.5-inch Hinged CryoBox with adjustable plastic grid, 81-place		2-inch Hinged CryoBox with adjustable plastic grid, LN ₂ drain holes, 81-place	
Qty per pack	5	50	5	50	6	30	5	50
○ white	BCS-206	BCS-207						
● blue	BCS-206B	BCS-207B						
● green	BCS-206G	BCS-207G	BCS-209G	BCS-220G	BCS-215G	BCS-219G	BCS-217G	BCS-221G
● orange	BCS-206O	BCS-207O						
● purple	BCS-206P	BCS-207P	BCS-209P	BCS-220P	BCS-215P	BCS-219P	BCS-217P	BCS-221P
● pink	BCS-206PK	BCS-207PK						
● multi-color (no white)	BCS-206MC							
Dimensions (L x W x H)	13.3 x 13.3 x 5.1 cm 5.25 x 5.25 x 2.0 in		13.3 x 13.3 x 5.1 cm 5.25 x 5.25 x 2.0 in		13.3 x 13.3 x 9.2 cm 5.25 x 5.25 x 3.5 in		13.3 x 13.3 x 5.1 cm 5.25 x 5.25 x 2.0 in	